

# MST 2007

## Measurement protocol



SALZMANN · GROUP  
ST. GALLEN · SWITZERLAND



**Group :** 醫碩科技股份有限公司


**Article :** 15-20mmHg Knee high socks

**Type :** A01

**Description :**

**Class :** I

**Size :** L/XL



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,9	19,2	1,4	1,6	0	0	100	96,5	7,0	8,0	0,0	0,0	02.08.2024 14:45:57
19,8	18,1	1,2	0,7	0	0	100	91,4	6,1	3,5	0,0	0,0	02.08.2024 14:52:33
19,9	18,7	1,3	1,2	0,0	0,0	100	94	6.5	5.8	0	0	Average

# MST 2007

## Measurement protocol



SALZMANN · GROUP  
ST. GALLEN · SWITZERLAND

**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A01

**Description :**

**Class :** I

**Size :** S/M



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
20	19,6	1,5	0,4	0	0	100	98,0	7,5	2,0	0,0	0,0	02.08.2024 15:41:00
19,5	19	0,6	1	0	0	100	97,4	3,1	5,1	0,0	0,0	02.08.2024 16:27:48
19,8	19,3	1,1	0,7	0,0	0,0	100	97,7	5,3	3,5	0	0	Average

# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A02

**Description :**

**Class :** I

**Size :** L/XL



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,9	19,6	0,8	0,9	0	0	100	98,5	4,0	4,5	0,0	0,0	02.08.2024 14:53:29
19,6	19,2	0,9	0,9	0	0	100	98,0	4,6	4,6	0,0	0,0	02.08.2024 16:28:46
19,8	19,4	0,9	0,9	0,0	0,0	100	98.2	4.3	4.6	0	0	Average

# MST 2007

## Measurement protocol



SALZMANN · GROUP  
ST. GALLEN · SWITZERLAND

**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A02

**Description :**

**Class :** I

**Size :** S/M



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,6	19,6	1,1	1,1	0	0	100	100,0	5,6	5,6	0,0	0,0	02.08.2024 15:43:35
19,5	18,8	0,9	0,8	0	0	100	96,4	4,6	4,1	0,0	0,0	02.08.2024 16:29:34
19,6	19,2	1,0	1,0	0,0	0,0	100	98.2	5.1	4.9	0	0	Average

## Measurement protocol



**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A03

**Description :**

**Class :** I

**Size :** L/XL



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
20	19,4	1,2	0,8	0	0	100	97,0	6,0	4,0	0,0	0,0	02.08.2024 14:57:12
19,1	19,1	0,9	0,7	0	0	100	100,0	4,7	3,7	0,0	0,0	02.08.2024 15:00:46
19,6	19,3	1,1	0,8	0,0	0,0	100	98.5	5.4	3.8	0	0	Average

# MST 2007

## Measurement protocol



SALZMANN · GROUP  
ST. GALLEN · SWITZERLAND



Group : 醫碩科技股份有限公司

Article : 15-20mmHg Knee high socks

Type : A03

Description :

Class : I

Size : S/M

Examiner :

Signature :



b b1 c d f g

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,8	19,3	1,1	0,9	0	0	100	97,5	5,6	4,5	0,0	0,0	02.08.2024 15:49:04
19,6	19	0,8	0,4	0	0	100	96,9	4,1	2,0	0,0	0,0	02.08.2024 15:52:58
19,7	19,2	1,0	0,7	0,0	0,0	100	97.2	4.8	3.3	0	0	Average

# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A04

**Description :**

**Class :** I

**Size :** L/XL



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,7	19,7	1,3	1	0	0	100	100,0	6,6	5,1	0,0	0,0	02.08.2024 15:01:44
19,6	19,6	1,1	0,8	0	0	100	100,0	5,6	4,1	0,0	0,0	02.08.2024 15:06:39
19,7	19,7	1,2	0,9	0,0	0,0	100	100	6.1	4.6	0	0	Average

# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司


**Article :** 15-20mmHg Knee high socks

**Type :** A04

**Description :**

**Class :** I

**Size :** S/M



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
20	19,5	1,1	0,4	0	0	100	97,5	5,5	2,0	0,0	0,0	02.08.2024 15:54:18
19,5	19,5	1,2	1,1	0	0	100	100,0	6,2	5,6	0,0	0,0	02.08.2024 15:55:00
19,8	19,5	1,2	0,8	0,0	0,0	100	98,7	5,8	3,8	0	0	Average



# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A05

**Description :**

**Class :** I

**Size :** L/XL



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
20	20	1	1,1	0	0	100	100,0	5,0	5,5	0,0	0,0	02.08.2024 15:08:13
19,9	19,9	1,1	0,8	0	0	100	100,0	5,5	4,0	0,0	0,0	02.08.2024 15:10:17
20,0	20,0	1,1	1,0	0,0	0,0	100	100	5.3	4.8	0	0	Average

Measurement protocol



Group : 醫碩科技股份有限公司

Article : 15-20mmHg Knee high socks

Type : A05

Description :

Class : I

Size : S/M



b b1 c d f g

Examiner :

Signature :

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,4	19,4	1,2	0,8	0	0	100	100,0	6,2	4,1	0,0	0,0	02.08.2024 15:56:08
19,2	18,5	1,5	1,1	0	0	100	96,4	7,8	5,7	0,0	0,0	02.08.2024 15:57:38
19,3	19,0	1,4	1,0	0,0	0,0	100	98.2	7	4.9	0	0	Average

# MST 2007

## Measurement protocol



SALZMANN · GROUP  
ST. GALLEN · SWITZERLAND

**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A06

**Description :**

**Class :** I

**Size :** L/XL



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,7	19,7	1	1,2	0	0	100	100,0	5,1	6,1	0,0	0,0	02.08.2024 15:14:41
20	20	0,8	0,5	0	0	100	100,0	4,0	2,5	0,0	0,0	02.08.2024 15:18:58
19,9	19,9	0,9	0,9	0,0	0,0	100	100	4.5	4.3	0	0	Average

# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A06

**Description :**

**Class :** I

**Size :** S/M



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19	19	1,3	1,3	0	0	100	100,0	6,8	6,8	0,0	0,0	02.08.2024 15:58:21
19,3	18,6	0,8	0,8	0	0	100	96,4	4,1	4,1	0,0	0,0	02.08.2024 15:59:07
19,2	18,8	1,1	1,1	0,0	0,0	100	98.2	5.5	5.5	0	0	Average

# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A07

**Description :**

**Class :** I

**Size :** L/XL



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
20,1	19,7	0,7	0,8	0	0	100	98,0	3,5	4,0	0,0	0,0	02.08.2024 15:20:05
19,9	19,9	1,1	0,9	0	0	100	100,0	5,5	4,5	0,0	0,0	02.08.2024 15:22:05
20,0	19,8	0,9	0,9	0,0	0,0	100	99	4.5	4.3	0	0	Average

# MST 2007



## Measurement protocol



SALZMANN · GROUP  
ST. GALLEN · SWITZERLAND

**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A07

**Description :**

**Class :** I

**Size :** S/M



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,3	18,6	0,9	0,8	0	0	100	96,4	4,7	4,1	0,0	0,0	02.08.2024 16:01:27
19,1	19,1	1	0,9	0	0	100	100,0	5,2	4,7	0,0	0,0	02.08.2024 16:02:13
19,2	18,9	1,0	0,9	0,0	0,0	100	98,2	4,9	4,4	0	0	Average

# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A08

**Description :**

**Class :** I

**Size :** L/XL



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,9	19,9	0,6	0,4	0	0	100	100,0	3,0	2,0	0,0	0,0	02.08.2024 15:23:56
19,9	19,9	0,8	1	0	0	100	100,0	4,0	5,0	0,0	0,0	02.08.2024 15:24:37
19,9	19,9	0,7	0,7	0,0	0,0	100	100	3.5	3.5	0	0	Average

# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A08

**Description :**

**Class :** I

**Size :** S/M

**Examiner :**

**Signature :**



b b1 c d f g

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,5	18,2	1,2	0,7	0	0	100	93,3	6,2	3,6	0,0	0,0	02.08.2024 16:03:32
19,4	18,6	1	0,7	0	0	100	95,9	5,2	3,6	0,0	0,0	02.08.2024 16:04:21
19,5	18,4	1,1	0,7	0,0	0,0	100	94,6	5,7	3,6	0	0	Average



# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司


**Article :** 15-20mmHg Knee high socks

**Type :** A09

**Description :**

**Class :** I

**Size :** L/XL



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
19,9	19,4	1,1	1	0	0	100	97,5	5,5	5,0	0,0	0,0	02.08.2024 15:28:17
20	19,4	0,4	0,8	0	0	100	97,0	2,0	4,0	0,0	0,0	02.08.2024 15:29:07
20,0	19,4	0,8	0,9	0,0	0,0	100	97,2	3,8	4,5	0	0	Average

# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A09

**Description :**

**Class :** I

**Size :** S/M



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
18,9	18,9	1,1	0,8	0	0	100	100,0	5,8	4,2	0,0	0,0	02.08.2024 16:05:05
19,6	18,6	1,1	0,7	0	0	100	94,9	5,6	3,6	0,0	0,0	02.08.2024 16:16:52
19,3	18,8	1,1	0,8	0,0	0,0	100	97.4	5.7	3.9	0	0	Average

# MST 2007

## Measurement protocol



SALZMANN · GROUP  
ST. GALLEN · SWITZERLAND

**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A10

**Description :**

**Class :** I

**Size :** L/XL



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
20	19,5	0,9	0,9	0	0	100	97,5	4,5	4,5	0,0	0,0	02.08.2024 15:33:03
19,6	19,6	0,6	0,6	0	0	100	100,0	3,1	3,1	0,0	0,0	02.08.2024 15:38:55
19,8	19,6	0,8	0,8	0,0	0,0	100	98.7	3.8	3.8	0	0	Average

# MST 2007

## Measurement protocol



**Group :** 醫碩科技股份有限公司

**Article :** 15-20mmHg Knee high socks

**Type :** A10

**Description :**

**Class :** I

**Size :** S/M



b b1 c d f g

**Examiner :**

**Signature :**

Pressure [mmHG]						Remaining Pressure [%]						
b	b1	c	d	f	g	b	b1	c	d	f	g	
20	19,5	1,3	0,9	0	0	100	97,5	6,5	4,5	0,0	0,0	02.08.2024 16:22:30
19,6	19,1	1,1	0,8	0	0	100	97,4	5,6	4,1	0,0	0,0	02.08.2024 16:26:48
19,8	19,3	1,2	0,9	0,0	0,0	100	97,5	6,1	4,3	0	0	Average